

Manual Supplement

Manual Title: BP Pump 2 Operator Supplement Issue: **2**
Part Number: 2196592 Issue Date: 9/10
Print Date: June 2007 Page Count: 1
Revision/Date:

This supplement contains information necessary to ensure the accuracy of the above manual.

FLUKE®

Biomedical

Change #1, 55935

On page 1-11, following **Pressure Generation** add:

Dynamic NIBP Simulation Repeatability..... ± 2 mmHg (at maximal pulse size independent of device under test)

Change #2, 56150, 56170

On page 1-11, under **Pressure Measurement**, change the High-Accuracy Version to the following:

High-Accuracy Version..... ± 0.7 mmHg (0.09 kPa)